

# Neuropixels Toolkit Live Script

## Build and execute a processing session with pipelines

### Setup Logger

```
import npxtoolkit.internal.thirdparty.logging.log4m
logger = log4m.getLogger("npx.log");
logger.clearLog();
logger.setLogLevel(logger.DEBUG);
```

### Processing Session

```
import npxtoolkit.session.Session
import npxtoolkit.pipeline.Pipeline
import npxtoolkit.internal.config.Config
% define session
session = Session('Session Info');
```

### Setup Python Env

```
PYENV_PATH = '/home/ubuntu/anaconda3/envs/npx/bin/python';
session.setPyEnv(PYENV_PATH);
% init session by pipelines, stages and jobs
```

### Method 1 - Auto-assembled pipeline

```
%% Auto-assembled pipeline
pipeline = Pipeline('Pipeline0 Info', "configs/test_config.json");
pipeline.autoAssemble();
session.addPipeline(pipeline);
session.parExecute();
```

### Method 2 - Manually-asmubled pipeline

```
import npxtoolkit.pipeline.Pipeline
import npxtoolkit.stage.Stage
import npxtoolkit.tasks.CatGT
import npxtoolkit.tasks.KiloSort
import npxtoolkit.tasks.TPrime

%% Pipeline 1
pipeline1 = Pipeline('Pipeline0 Info', "configs/test_config.json");

% CatGT stage
stageCatgt = Stage('CatGT');
pipeline1.addStage(stageCatgt);
% CatGT task 1
taskCatgt1 = CatGT('CatGT probe 0', '0', 1, pipeline1.PipelineConfigs);
stageCatgt.addTask(taskCatgt1);

% KiloSort stage
stageKilo = Stage('KiloSort');
```

```

pipeline1.addStage(stageKilo);
% KiloSort task 1
taskKilo1 = KiloSort('KiloSort probe 0', '0', 'cortex', pipeline1.PipelineConfigs);
stageKilo.addTask(taskKilo1);

% TPrime stage
stageTPrime = Stage('TPrime');
pipeline1.addStage(stageTPrime);
% TPrime task 1
taskTPrime1 = TPrime('TPrime probe 0', '0', pipeline1.PipelineConfigs);
stageTPrime.addTask(taskTPrime1);

% append pipeline to session
session.addPipeline(pipeline1);
%% Execution
session.parExecute();

```

## Single Task Execution for testing or debugging

```

import npxtoolkit.internal.config.PipelineConfig
import npxtoolkit.tasks.CatGT
import npxtoolkit.tasks.KiloSort
import npxtoolkit.tasks.TPrime

pipelineConfig = PipelineConfig("configs/test_config.json");
task = CatGT('CatGT probe 0', '0', 1, pipelineConfig);
% config = TaskConfig(json.KiloSort);
% task = KiloSort('KiloSort probe 0', '0', 'cortex', pipelineConfig);
% config = TaskConfig(json.TPrime);
% task = TPrime('TPrime probe 0', '0', pipelineConfig);
task.execute();

```

## Visualization Examples

These are currently not included in the Neuropixels-toolkit package, just to demonstrate the pipeline's results.

Data files are coming from the data directory in the config json file config["data"], in this case, they are in: `/home/ubuntu/neuropixel/data\_for\_ecephys/subject1\_session1/catgt\_SC011\_021919\_g0/SC011\_021919\_g0\_imec0/imec0\_ks2`

```

% this is the same as config["pipeline"]["tools"]["numpyMatlabRepo"]
addpath('/home/ubuntu/neuropixel/np-matlab/np-matlab/');
addpath('.');

```

## Mean Waveforms

```

meanWaveforms = readNPY('./test_data/mean_waveforms.npy');
disp(size(meanWaveforms));

```

```

275    301    82

```

```

samplingRate = 29999.9

```

```


samplingRate = 3.0000e+04

```

```

one_unit = squeeze(meanWaveforms(2,:,:))';
plot((0:size(one_unit, 1)-1)/samplingRate*1000, one_unit);
title('Mean Waveforms');
xlabel('Time(ms)');
ylabel('Amp');

```



## Spike Time

```

spikeTime = readNPY('./test_data/spike_times.npy');
disp(size(spikeTime));

```

```

7971090      1

```

```

spikeCluster = readNPY('./test_data/spike_clusters.npy');
disp(size(spikeCluster));

```

```

7971090      1

```

```

scatter(spikeTime(1:10000:7971090), spikeCluster(1:10000:7971090));
title('Spike Time vs Spike Clusters downsampled by 10000');
xlabel('Spike Time');
ylabel('Spike Cluster');

```

